## **Quadrupole Miniaturization – Reconsidered**

Justin Jacobsky<sup>1</sup>, Randall E. Pedder<sup>2</sup>, Michael Klaric<sup>2</sup>, Christopher Taormina<sup>2</sup>,

<sup>1</sup>Department of Chemical and Petroleum Engineering, University of Pittsburgh, Pittsburgh, PA; <sup>2</sup>Ardara Technologies L.P., Ardara, PA

Eight years ago, at the 2007 HEMS workshop, we presented our work in characterizing quadrupole performance as quadrupoles are miniaturized. We concluded that presentation with the statement: "...extreme miniaturization of quadrupoles is unlikely to yield analyzers with anywhere near the performance of conventionally scaled quadrupoles." Maybe we were a little hasty...

This year we re-visited quadrupole performance for not-so-extreme miniaturization, using quadrupoles with reduced rod length (4 cm), and reduced rod diameter (4 mm), operating at 6.5 MHz, with a mass range of 100 amu, with unexpectedly remarkable results. We developed a miniaturized closed ion source, scaled to optimize coupling with this small quadrupole, with a glass capillary inlet to sample atmospheric air or other samples brought to a simple atmospheric pressure Flow-by-Tee inlet, and adapted an off-the-shelf 2 <sup>3</sup>/<sub>4</sub> inch Faraday/Multiplier detection flange (from Extorr), all mounted to a five-inch-long vacuum chamber, pumped by a Pfeiffer Vacuum Hi-Pace 10 pump. Reducing quadrupole length promises to allow operation at higher than normal operating pressures; this work evaluates performance at pressures as high as one millitorr chamber pressure.

First results from this characterization show better than one ppm limit-of-detection for a single five-second-long, 50 amu profile spectrum (0-50 amu, 10 points-per-amu, oversampled with a 100 KHz data acquisition rate). Helium in room air is clearly visible at m/z 4, very near the limit-of-detection for this analyzer under these conditions. Traditional mass analyzer figures of merit (abundance sensitivity, limit-of-detection, mass resolution...) will be presented for this quadrupole system, along with representative spectra at various pressures, allowing evaluation of the limitations of this small analyzer.