

OPTIMISATION OF A PORTABLE QMS FOR ENVIRONMENTAL MONITORING

T. Fildes, D.G. McIntosh, M.J. Antony Joseph, B. Smith, S. Maher,
J.R. Gibson and S. Taylor

Mass Spectrometry Group
Department of Electrical Engineering and Electronics,
University of Liverpool, U.K.



UNIVERSITY OF
LIVERPOOL


TRACKING ENVIRONMENTAL CHANGE

- Global need for environmental data
- Methods of capturing environmental information *in situ* are few, and limited
- Field analysis restricted by technique
- High reliance on laboratory for complex analysis
- Geographic spread of datasets restricted by technology (cost and portability)



TRACKING ENVIRONMENTAL CHANGE

- Stable isotope measurements (IRMS) inform global trends
- e.g. $\delta^{13}\text{C}$ falls with increasing CO_2
- Small fluctuations significant
- $$\delta(\text{‰}) = \left(\frac{R_{\text{sample}} - R_{\text{standard}}}{R_{\text{standard}}} \right) \times 1000$$
- High measurement precision required: better transmission \rightarrow greater precision
- Other applications: archaeology < 2 (‰) paleobiology, volcanology, geochronology...



THE NEED FOR PERFORMANCE

Mass Spectrometry

- FT-ICR

MALDI-TOF


- Orbitrap
- AMS



THE NEED FOR COMPACT

Quadrupole Mass Spectrometry

- Versatile, scalable and accessible
- Simple design and affordable
- BUT Smaller → worse performance (usually)




- Can we improve small footprint QMS performance to allow on-site IRMS?

QMS OPTIMISATION

Variables...

- Filter geometry
 - Rod length
 - r_0
 - r/r_0
 - Source gap
 - r_{ie}
- Pre/post-filter(s)
- DC offset
- Rod bias



Goals...

- Efficient **transmission**
- Optimal **stability** (precision)

whilst maintaining

- Acceptable **resolution**
- Compact design
- Affordability


METHODS


Simulations in x, y, z and time supported by experimental measurements

Source Gap – space between source and QMF (or pre-filter)

Pre-filter – RF-only to aid transmission efficiency /stability


Rod Bias – DC voltage applied to QMF (and/or pre-filter) electrodes is raised to slow down ions from the source





SIMULATION TOOL QMS-3D: Stable ion trajectories 3D animation

Screen capture from 'QMS_Traj' (written in MATLAB: D. McIntosh, Mass Spectrometry Group, UoL, UK)




Unstable trajectory 3D animation demo

Screen capture from 'QMS_Traj' (*written in MATLAB: D. McIntosh, Mass Spectrometry Group, UoL, UK*)

1. SOURCE GAP




- 3D fringe fields at QMF entrance affect ion transmission: ion motion coupled in x, y and z
- Prior experimental treatments: fringe fields reduce ion transmission efficiency for heavier ions [1,2]
- Heavier ions spend a longer period exposed to defocusing forces in the y direction [1,2]



1. Brubaker WM. Auxiliary electrodes for quadrupole mass filters. US Patent 3129327. 14 Apr. 1964
2. Brubaker WM. An improved quadrupole mass analyzer. *Adv. Mass Spectrom* 1968; 4293-299.

2. SOURCE GAP


- **Ehlert**: isotope ratios → experimental **transmission curve** [1]
- Loss at low m/z ... “**artefact**” [1]? BUT our experiments (single filter QMS) show both rise and fall are real [2]
- Loss at lower m/z is due to **fringe field effects**
- New method to calculate **absolute experimental transmission**
- QMS-3D simulation model [2, 3] accurately predicts optimum transmission location (in m/z) and percentage



1. Ehlert TC. Determination of transmission characteristics in mass filters. *J.Phys E: Scientific Instr.* 1970; 3(3):237-239.
2. Antony Joseph M, McIntosh D, Gibson R, Taylor S. Effects of the source gap on transmission efficiency of a quadrupole mass spectrometer. *RCMS*, 2018.
3. Gibson JR, Evans KG, Syed SU, Maher S, Taylor S. A method of computing accurate 3D fields of a quadrupole mass filter and their use for prediction of filter behavior. *J.American Society for Mass Spectrometry* 2012; 23(9):1593-1601.



2. SOURCE GAP

- Length of source gap dictates fringe field length and field intensity
- Not analogous to altering ion energy
- QMS-3D model allows direct study of relationship between physical source gap length and stable ion transmission efficiency
- Different source gap lengths promote peak transmission in different m/z ranges [1]
- **Source gap can be tailored for optimum transmission in desired mass range**



1. Antony Joseph M, McIntosh D, Gibson R, Taylor S. Effects of the source gap on transmission efficiency of a quadrupole mass spectrometer. *Rapid Communications in Mass Spectrometry* 2018.


2. SOURCE GAP



- Smaller source gaps \rightarrow optimal transmission at higher m/z and vice versa [1]
- E.g 1.5r0 source gap improves transmission at m/z 12, but for m/z 120, 0.25 r0 source gap is better
 - Transmission gains are not at the expense of resolution
- Source gap optimisation: a critical component of QMS design for mass-specific applications [2,3]
 - **Accurate simulation of the specific instrument is key**

1. Antony Joseph M, McIntosh D, Gibson R, Taylor S. Effects of the source gap on transmission efficiency of a quadrupole mass spectrometer. *RCMS* 2018.
2. Ellefson RE, Moddeman WE, Dylla HF. Hydrogen Isotope Analysis by Quadrupole Mass-Spectrometry. *Journal of Vacuum Science & Technology* 1981; 18(3):1062-1067.
3. Schneider B, Kuiper K, Postma O, Wijbrans J. Ar-40/Ar-39 geochronology using a quadrupole mass spectrometer. *Quaternary Geochronology* 2009; 4(6):508-516.


UPDATE: BRIEF PROGRESS IN $\delta^{13}\text{C}$ STABILITY






- Example of progress excluding current optimisation study
- **Dual Inlet:** well established standard stable isotope apparatus
 - Subtraction of temporal instability

$$\delta(\text{‰}) = \left(\frac{R_{\text{sample}} - R_{\text{standard}}}{R_{\text{standard}}} \right) \times 1000$$

- *Instruments sold / in service*
www.q-technologies.co.uk




PRE-FILTER EFFECT?



- With crimped capillary inlet, sample CO₂ @ ~34 psi
- Comparison not strictly fair: different geometries modes /settings of acquisition/ on-board signal processing...
- Could **pre-filter** be a factor? Simulation a fairer test

2. PRE-FILTER

- Brubaker used a delayed D.C. ramp (pre-filter) to reduce the effects of source gap fringe fields [1, 2]
- Improved transmission of low velocity ions (high m/z)
- Using QMS-3D simulation, can design a minimum prefilter length for given electrode, frequency and ion energy [3]
- Caution: for low m/z ions (high velocity), see transmission 'fluctuations' (function of ion velocity)



Ion transmission for several lengths of prefilter


References

1. Brubaker WM. Auxiliary electrodes for quadrupole mass filters. *US Patent* 3129327. 14 Apr. 1964
2. Brubaker WM. An improved quadrupole mass analyzer. *Adv. Mass Spectrom* 1968; 4293-299.
3. Gibson JR, Evans KG, Taylor S. Predicted behaviour of QMF systems with and without prefilters using accurate 3D fields. *Int J. of Mass Spec.* 2017; 422197-207.

QMS-3D Simulations

- With and without prefilter
- 20 scans each, 20,000 ions/scan
- Upper series: with prefilter
- Transmission: 23-24% improvement
- *Slight* improvement in resolution (FW @ 10%)
- Effect on quantitative measurement stability?

2. PRE-FILTER



Simulations: stability predictions

	Single filter	Dual filter	Effect
Average Trans.	21.8%	26.9%	23% better
Standard deviation	0.23	0.17	24% better
RSD	1.03%	0.64%	38% better

Tracking maximum transmission

- Pre-filter boosts transmission and improves quantitative measurement reproducibility
 - Worse standard deviation in the fixed m/z case could imply an effect on mass filtering accuracy ('fluctuations' relevance?? ... more study required)



2. PRE-FILTER

	Single filter	Dual filter	Effect
Average Trans.	21.57%	26.71%	24% better
Standard deviation	0.24	0.26	9% worse
RSD	1.12%	0.98%	12% better

Transmission at a fixed m/z point

3. ROD BIAS



- DC bias of all electrodes in one direction
 - Independent of rod initial mean polarity
- Rod bias alters mean potential gradient in entrance fringe field
- Can be applied:
 - directly to electrodes
 - inversely to entrance / exit plates
- Slower ions \rightarrow more RF cycles \rightarrow greater resolution for the same transmission (and vice versa) ?



Stability diagram - Dawson, Peter H., ed. *Quadrupole mass spectrometry and its applications*. Elsevier, 2013.

QMF electrical diagram - Brubaker, Wilson M. "Auxiliary electrodes for quadrupole mass filters." U.S. Patent No. 3,129,327. 14 Apr. 1964.

3. ROD BIAS









- Use easily variable SIMION ion optics to simulate rod bias effects
- Transmission trend plotted against “Axial Energy (Ae)” for varying rod bias
- At lower Ae, transmission with biased rods surpasses unbiased rods
- Effect is more pronounced at higher resolution (U/V ratios)

Low resolution study

- Rod bias \rightarrow low axial energies
- Rod bias improves transmission at extremely low axial energies
- Improvement in peak shape with added bias
- Resolution constant







3. ROD BIAS



Higher resolution study

3. ROD BIAS





- Rod bias: best at low axial energies
- Allows greater transmission without significantly compromising resolution



Experimental (single filter QMS)



- Axial energies by varying injection energies (ie) and rod bias values (rb)
- At lower Ae some rb:ie values allow improved transmission without resolution cost

3. ROD BIAS



Experimental v SIMION

3. ROD BIAS



- Ion injection energy of 10eV, increasing rod bias from 0-10V
- Experiments: steady incremental drop; simulations: more dramatic decline, as **rb** --> **ie**
 - Different parameters (e.g. source gap)? Unknown factors?

SUMMARY

Isotope Ratio MS (IRMS) needs high ion transmission to achieve stability required

1. Source gap

- Favours ion transmission efficiency in the (low) m/z range for $\delta^{13}C$ and $\delta^{15}N$
- No cost to resolution

2. Prefilter

- Improves ion transmission efficiency and measurement stability
- No cost to resolution

3. Rod bias

- Certain combinations of rod bias and ion injection energy can improve transmission without compromising resolution

CONCLUSIONS AND FURTHER WORK



1-500 Da portable QMS

- Each method independently yields gains
 - Relevance of the post-filter?
- Challenge: thread the approaches together
 - Optimal design of a small footprint QMS
 - High-fidelity simulation is key to this goal
- Combine with dual inlet miniaturisation for a portable IR-QMS system for onsite IRMS for
 - Environmental monitoring (VOCs), volcanology, archaeology, paleobiology, and geochronology

IMAGE CREDITS

Oil burning - <http://mansfieldgaswellawareness.weebly.com/industry-news.html>

Worker taking readings - <https://www.ventia.com/>

Laboratory IRMS - <https://www.marine.usf.edu/research/facilities-and-equipment/stable-isotope-ratio-mass-spectrometry-facility/>

FT-ICR - <https://www.gfz-potsdam.de/en/section/organic-geochemistry/infrastructure/ft-icr-ms/>

MALDI-TOF – <https://www.bruker.com/products/mass-spectrometry-and-separations/maldi-toftof.html>

Orbitrap - <https://www.thermofisher.com/order/catalog/product/IQLAAEGAAPFALGMAZR>

AMS - <http://www.pelletron.com/products/accelerator-mass-spectrometry-ams-systems/>

QMF - <https://www.henniker-scientific.com/products/instruments/gas-phase-chemistry/high-resolution-quadrupole-mass-spectrometers>

[all other illustrations: authors]

ACKNOWLEDGEMENTS

We would like to acknowledge the contributions and advice of:


- Dr. K.G. Evans (University of Liverpool)
- Peter Chaisty and Milos Dragovic (Cyionics, UK)
- Dr. Neil France (Q-Technologies, UK)



UNIVERSITY OF
LIVERPOOL

The background of the slide is a solid black field. At the top, there is a decorative horizontal band with a wavy, fluid appearance. This band features a color gradient: on the left, it transitions from yellow to orange; on the right, it transitions from green to cyan. The colors blend into each other, creating a sense of motion and depth.

THANK YOU! QUESTIONS?



EXTRA SLIDES IF REQUIRED


- 1st: for background / reminder of HySEM-80 early raw d13C results without prefilter, without crimped capillary, etc. (From HEMS last yr)
- 2nd: for some old VapourSense results on N2 isotope stability – if required to satisfy the abstract. (These were presented at ASMS 2017 and HEMS last yr)
- Both could go before slide 11
- PS if preferred,: Slide 11 could be moved to be placed after slide 3 before getting into the main approaches in today's presentation

STABILITY AND LINEARITY – EARLY RESULTS: SINGLE FILTER HYSEM-80

CO₂ gas various pressures 1



CO₂ gas various pressures 1



<input checked="" type="checkbox"/>	loR_1.89A_7.9e-6_hIQ_CO2_gas (2)
<input checked="" type="checkbox"/>	loR_1.89A_7.9e-6_hIQ_CO2_gas (3)
<input checked="" type="checkbox"/>	loR_1.89A_7.9e-6_hIQ_CO2_gas (4)
<input checked="" type="checkbox"/>	loR_1.89A_7.9e-6_hIQ_CO2_gas (5)
<input checked="" type="checkbox"/>	loR_1.89A_4.0e-6_hIQ_CO2_gas
<input checked="" type="checkbox"/>	loR_1.89A_3.5e-6_hIQ_CO2_gas (2)
<input checked="" type="checkbox"/>	loR_1.89A_3.5e-6_hIQ_CO2_gas (3)
<input checked="" type="checkbox"/>	loR_1.89A_3.5e-6_hIQ_CO2_gas (4)
<input checked="" type="checkbox"/>	loR_1.89A_3.5e-6_hIQ_CO2_gas (5)
<input checked="" type="checkbox"/>	loR_1.89A_3.2e-6_hIQ_CO2_gas
<input checked="" type="checkbox"/>	loR_1.89A_3.2e-6_hIQ_CO2_gas (1)
<input checked="" type="checkbox"/>	loR_1.89A_2.8e-6_hIQ_CO2_gas
<input checked="" type="checkbox"/>	loR_1.89A_2.8e-6_hIQ_CO2_gas (1)
<input checked="" type="checkbox"/>	loR_1.89A_2.5e-6_hIQ_CO2_gas
<input checked="" type="checkbox"/>	loR_1.89A_2.5e-6_hIQ_CO2_gas (1)
<input checked="" type="checkbox"/>	loR_1.89A_2.5e-6_hIQ_CO2_gas (1)

Measurements from line 1 of Table 1. Six clear groups can be discerned in the mass 44 peak (= mostly ¹²CO₂) corresponding to 6 reducing pressure bands of introduction of cylinder CO₂ gas to the needle valve inlet of the portable QMS. The first four of these can be clearly discerned on the inset (mass 45 = ~93% ¹³CO₂).

CO₂ cylinder gas

Pressure range (Torr)

Mass 45/44
standard
deviation

Predicted δ¹³C
standard
deviation

Linearity across 6 pressures, ~5 per
pressure

7.9 x 10⁻⁶ – 1.6 x 10⁻⁶

0.000252

22.4‰

Above – averaging groups of ~5 data
points

7.9 x 10⁻⁶ – 1.6 x 10⁻⁶

0.000159

14.1‰

Above – averaging groups of ~5 ratios

7.9 x 10⁻⁶ – 1.6 x 10⁻⁶

0.000122

10.9‰

Above – with high outlying pressure
removed

3.6 x 10⁻⁶ – 1.6 x 10⁻⁶

0.000089

7.84‰

Above – with high outlying pressure
removed

3.6 x 10⁻⁶ – 1.6 x 10⁻⁶

0.000090

7.99‰


- Need turbulent flow at inlet
- Standard to use a capillary
- Reynold's number: $Re = \frac{\bar{v}D\rho}{\mu}$
 \bar{v} = mean velocity
 ρ = gas density
 D = vessel diameter
 μ = gas viscosity

- ***Re* needs to exceed critical value > ~2000 (ideally higher) ... 0.15mm I.D. capillary → crimped**

STABILITY – PRIOR VAPOURSENSE-500 RESULTS (N₂)

- Peak-jumping mode (MIM)
- Sample not pressurised
- Capillary not standard
- Sintered leak may cause fractionation
- 29/28 signal is surprisingly stable despite non-ideal setup

VapourSense-500, ambient N ₂ MIM mode (peak jumping)	m/z 29/28 stand. dev. (1 σ)
Stability across 247 measurements (34-minute period)	0.0000501
Above – final 37 measurements only (5- minute period)	0.0000275



METHODS

Simulations in x , y , z
and time supported by
experimental
measurements

