

Two-Dimensional Tandem Mass Spectrometry for Bacterial Profiling


L. Edwin Gonzalez, Lucas J. Szalwinski, Thomas C. Sams, Eric T. Dziekonski, R. Graham Cooks

Department of Chemistry, Purdue University
Aston Laboratory for Mass Spectrometry Research



Introduction

- Motivation
- Sample Preparation
- Ion Trap Theory
- 2D MS/MS Analysis
- 2D Opportunities



Biological Threats





Table 1: Agents that can be used as biological warfare or bioterrorism-related incidents

	Agents	Disease	Route of infection	Possible release
Bacteria	<i>Ba. anthracis</i>	Anthrax	Aerosol	Spores
	<i>Y. pestis</i>	Plague	Aerosol	Vegetative cells
	<i>Br. melitensis</i>	Brucellosis	Aerosol	Vegetative cells
	<i>Br. abortus</i>			
	<i>Bu. mallei</i>	Glanders	Aerosol	Vegetative cells
	<i>Bu. pseudomallei</i>	Melioidosis	Aerosol	Vegetative cells
Viruses	Variola virus	Smallpox	Aerosol	Virus particles
	Ebola virus	Ebola hemorrhagic fever	Aerosol	Virus particles
Toxins	Marburg virus	Marburg hemorrhagic fever	Aerosol	Virus particles
	<i>C. botulinum</i>	Botulism	Ingestion food/water	Toxin
	<i>Staphylococcus aureus</i>	Staphylococcal enterotoxin type B (SEB)	Food/water	Toxin
	Ricin (plant)	Ricin toxin	Food/water	Toxin
	Trichothecene (fungus)	Trichothecene T2 toxin	Food/water	Toxin

1) <https://www.cdc.gov/anthrax/basics/types/index.html#:~:text=Inhalation%20anthrax%20starts%20primarily%20in,about%2055%25%20of%20patients%20survive>.

2) Thavaselvam D, Vijayaraghavan R. Biological warfare agents. *J Pharm Bioallied Sci.* 2010 Jul;2(3):179-88. doi: 10.4103/0975-7406.68499. PMID: 21829313; PMCID: PMC3148622.

3) https://en.wikipedia.org/wiki/2001_anthrax_attacks

Bacterial Profiling


- *Bacillus* Genus Form Spores

- Each Species has a Unique Lipid Membrane


- Dipicolinic Acid is a Metabolic Product of Sporulation

- Lyse Spores with Microwave

Phosphatidylglycerol (PG)




Dipicolinic acid





Mass Spectrometry


- nESI
- Modified Thermo LTQ LX
- Supplementary Waveforms
- Linear Quadrupole Ion Trap
- MSSuite Software





Two-Dimensional Tandem Mass Spectrometry






Two-Dimensional Tandem Mass Spectrometry





Sporulation Metabolite of *B. subtilis*





Sporulated Lipids of *B. subtilis*





Vegetative Metabolite of *B. subtilis*





Vegetative Lipids of *B. subtilis*



Sporulated Lipids of *B. thurgiensis*



Bacterial Identification



Opportunities

- Threat Analysis
- Lipidomics
- Lipid Libraries
- Rapid In Situ Analysis
- Miniaturization
- Qualitative Analysis
- Harsh Environments



Acknowledgments

- Thanks to Edwin, Lucas, Eric, and Graham
- Thanks to the rest of Aston Labs
- Funded by the Department of Homeland Security

