

Chip Scale Mass Spectrometry

...Yet another contender for Harsh Environments? Or can it help redefine "Harsh Environments"?

Presented at – The 14th Harsh-Environment Mass Spectrometry Workshop
September 28, 2022

Ashish Chaudhary (Ph.D.)

(Founder and CEO)

Detect-Ion LLC

Ashish.chaudhary@detect-ion.com

Tel: (727) 251-0889



3802 Spectrum Blvd. (Suite 128)


Tampa, FL 33612

www.detect-ion.com

This research is based on work supported by the Office of the Director of National Intelligence (ODNI), Intelligence Advanced Research Projects Activity (IARPA), via FA8650-17-C-9108, FA8650-19-C-9100 and 20120800003. The views and conclusions contained herein are those of the authors and should not be interpreted as necessarily representing the official policies or endorsements, either express or implied, of the ODNI, IARPA, or the U.S. Government. The U.S. Government is authorized to reproduce and distribute reprints for Governmental purposes notwithstanding any copyright annotation thereon.

Acknowledgments

Ryan Bell, Friso van Amerom, John Nicosia, Judd Sheets, James Baughman, Wei yang and more!



SRI International®

Matt Colvin, Strawn Toler, Tim Short, Dustin McRae, Jim Alberti, John Kelly, Jeff Strickrott, Michael Isnardi, Michael Wei, Billy Hoffmann, Annie Dowgiallo, and Ashish Chaudhary



Detect-Ion LLC

("Detection")

Mission: Next Generation Sensor Development For Enhanced CBRNE Sensing

Vision: Bring Traditional Mass Spectrometry Performance to the Masses


- Small business technology startup (October 2021)
- Technical team hails from SRI International
 - Performers on IARPA MAEGLIN and DARPA SIGMA+ CWMD programs
- Inducted into the Tampa Bay Technology Incubator (TBTI) and housed on USF Tampa campus
 - Dedicated laboratory space
 - Onsite rapid prototyping
 - Access to shared high-value lab resources through the TBTI program
 - Florida Hi-Tech Corridor Matching Funds




Proud partners with



Motivation: Disrupting the Monotonic SWaP vs. Performance Trade off



Research Themes that enabled ACHILLES






Roadmap

Half a lab!

Lots of ancillaries

17.4 L

~10 L



Validation Testbed (2017)



Lab Model (2018)



Transportable Prototype (2019)



ACHILLES (2021)


MAEGLIN P-1

SIGMA+ (Yr-1)

MAEGLIN P-2



ACHILLES Block Diagram



ACHILLES Prototype

Currently outside the ACHILLES box

LabVIEW-based control software


ACHILLES (TRL-5/6)

Technology Highlights

- World's lowest Swap-C [TD-GC-MS] platform
- High sensitivity validated for >200 VOCs
- Chemical portfolio: semi-volatile and volatile organic compounds (VOCs) including, but not limited to, common solvents, precursors, TICs, TIMs, CWA surrogates, and explosives related compounds (ERCs).
- Detection algorithm leverages separation (RT) and near 1-amu resolution(m/z) mass spectral signature
- Similarity with NIST mass spectral libraries

Independent Government Validation

- 100% (+ true known) in NRL test campaign*
- Successfully flight-testing for USAF/AFRL
- Grand Central Station (NYC) under DHS CBT program






*DeWitt K. "Advancements in compact gas collection and analysis from IARPA's MAEGLIN program." CBRNE Sensing XXI. Vol. 11416. International Society for Optics and Photonics, 2020.

Ion Trap: Background

Mass-selective instability mode operation of RF ion traps


$$V_{rf} = \frac{mq_z r_0^2 \Omega^2}{4A_2 e}$$




Miniaturization of ion traps in micromachined high-density arrays has led to low-SWaP MS design

We achieved near 1-amu resolution

Good mass spectral resolution is key to accurate chemical ID in real-world clutter environment




~200 traps operated in unison; while maintaining 1-amu mass resolution for the collective signal from the array



200-pixel chip-scale ion trap array


High-density chip-scale ion trap array with distributed localized ion trapping across 200 traps overcomes the peak broadening due to space-charge and ion-neutral collision effects that results in poor resolution of high-pressure miniature ion traps.

Partial mass spectrum of Trichloroethylene



Chip Scale MS: Limits of Detection (LOD) (nanogram)

- Chip-Scale MS has been demonstrated for ~200 chemicals with diverse chemical classes
- These LODs are 10-15x better than currently fielded systems (typically 1-10 ng)



ACHILLES (ppb may be more interesting than ng)

	Sampling time (s)			LODs
	120 s	60 s	30 s	
Toluene	3.23	6.45	12.90	$\mu\text{g}/\text{m}^3$
	0.86	1.71	3.42	ppbv
DMMP	17.25	34.5	69	$\mu\text{g}/\text{m}^3$
	3.39	6.79	13.58	ppbv
DEMP	81.45	162.9	325.8	$\mu\text{g}/\text{m}^3$
	13.07	26.15	52.29	ppbv

- LODs of chemicals listed above are a good representation of the LOD range for the volatile organic compounds (such as toluene) and semi-volatile organic compounds (such as DMMP and DEMP).
- ACHILLES can be operated with a wide range of sampling time to fit with the detection requirements.
- Concentrations higher than the range mentioned above will be detectable with a slight compromise in the mass spectral resolution. Additionally, the sampling time can be optimized for the target range of concentration, as needed.

Our GC split throws away >96% collected molecules; Only 3-4% reach the MS for detection




Sneak peak: We have figured out how not to do that!

Library Development

Red: NIST mass spectra

Black: μ MS mass spectra

- μ MS produces mass spectra similar to NIST libraries
 - Rapid library development for >100s of chemicals
 - “Gotcha’s” due to minor deviations
 - In MAEGLIN we developed custom library/chemical
- ML model could enable an [“infinite library”](#) platform





ACHILLES: System stability → Quantification and Identification

Back-to-back experiment (n=9)
from a stable test environment



Average RSD (intensity) across all six chemical components ~ 10%

Good RT stability (Std dev < 0.25 s)



RT (s)	23.61	37.85	45.11	59.35	72.88	76.08
Run-1	14.40	4.93	4.37	21.00	16.90	3.75
Run-2	11.50	4.67	4.18	22.80	18.00	3.79
Run-3	10.60	4.77	4.06	24.00	18.40	3.71
Run-4	10.70	4.89	3.91	19.70	16.90	3.10
Run-5	12.40	4.46	3.65	18.50	16.90	3.95
Run-6	13.90	5.09	5.35	22.00	19.10	4.85
Run-7	12.30	5.09	4.52	23.60	17.80	4.78
Run-8	16.40	4.84	4.90	20.80	19.20	4.42
Run-9	12.70	5.38	3.87	19.60	16.50	4.01
STD	1.87	0.27	0.54	1.90	1.01	0.56
Avg	12.77	4.90	4.31	21.33	17.74	4.04
RSD	14.67	5.45	12.56	8.92	5.69	13.83

Average RSD (Intensity) across all six chemical components ~10%



Important for absolute and relative quantification of VOCs

Important for identification

Identification Algorithm

Approach/Status:

- 2D GC-MS data
- Pre-developed library cards
 - Curated identification based on “maximally disjoint” AMU peak positions
 - Some adducts (evolving)


Process

- Background subtraction
- Noise reduction (Savitzky-Golay, FFT etc.)
- Peak finding, centroiding
- RT and AMU window threshold optimization



Key achievements so far:

- Repeatable detection of co-eluting chemicals
- <1-min analysis time (~133 cards)
- Low P_{fa} (1 - 2%)

Full 2D Map of processed GC-MS data



Library cards examined w.r.t. RT




Deployment/Use-Case



Outdoor Sampling (NRL)



Tropical Bay (NRL)





Mobile Deployment (Boston)




Grand Central Station (NYC)



USAF/AFRL Flight Testing



Breath Diagnostics



Dehumidification stage

Water sample

ACHILLES as Water Sensor

USAF CBRN Flight Testing (1/3)





Office of the Director of National Intelligence

I A R P A

BE THE FUTURE




- Objective: Real time measurement of Methyl Salicylate (MeS) with different flight maneuvers to understand the concentration change in a low altitude flight mode
- ACHILLES was operated during the flight (including take-off and landing) and was co-located with COTS chemical sensors
- During test, the forward cockpit hatch and both paratrooper doors were opened to purge the plane of the test compound




USAF CBRN Flight Testing (2/3)

Methyl Salicylate (MeS) is a common chemical agent surrogate


NIST MeS Mass Spectrum



ACHILLES MeS Mass Spectrum




MeS Concentration Estimates vs. Local Time




USAF CBRN Flight Testing (3/3)

Major achievements/takeaways

- ACHILLES measured MeS concentration trends which corroborated with the flight maneuvers and matched with the trend measured by the co-located COTS PID (PPBRAE4 detector)
- During the tail end of the exercise (“soak”), ACHILLES was able to detect the onset of rise in concentration (MeS outgassing from aircraft surfaces due to tarmac temperature), which were missed by the COTS
- As mentioned by the activity lead, AFRL had flow COTS mass spectrometers in the past, but those platforms suffered from failure modes during the flight



Harsh Environments (Requirements)



- Low system latency requires intelligent inputs for optimal sensor operation (Avoid saturation but maintain sensitivity)
 - Can complementary ionization enable a shorter GC method?
- Sampling to final answer (deconvolution, identification, quantification etc.)
 - Can we make the cleaning cycle intelligent based on signal detected?
- Ability to expand the libraries and move to a “training free” paradigm
 - Can an AI/ML model bridge the gap b/w ACHILLES and NIST libraries?
- Long endurance strategic emplacements limits serviceability
 - We are investigating N_2 and eventually “ambient air” as the consumable

Summary

- ✓ Chip Scale Mass Spectrometry shows a good promise to improve the state-of-the-art of fieldable trace chemical sensing
- ✓ There is room for further miniaturization and ruggedization
 - ACHILLES is the result of first prototype integration
 - Vacuum pump solution is an under-explored research theme
- ✓ Beyond hardware optimizations, software and application-specific AI/ML algorithm will be key areas of development necessary to reach a minimal viable product
- ✓ Detect-Ion is diversifying ACHILLES in areas of water security, health diagnostics, environmental monitoring

Thank you!