

DEVELOPMENT OF A MINIATURE RECTILINEAR ION TRAP ARRAY WITH INDEPENDENTLY CONTROLLED CHANNELS


Paul Hendricks,¹ Jeff D. Maas,² William J. Chappell,² Rob Noll,¹ Zheng Ouyang,³
R. Graham Cooks¹

¹Dept. of Chemistry, Purdue University, IN;


²Dept. of Electrical Engineering, Purdue University, IN;

³Dept. of Biomedical Engineering, Purdue University, IN

Mass analyzer miniaturization: Scale, RF, Fabrication Technologies, and Cost





- Complementary to miniaturization
- Fabrication of devices with complex geometries
 $y_0 = 1.33\text{mm}$
- Cost effective (\$1500 vs. \$110)
- Manufacturing time (days vs. weeks)
 $x_0 = 1.66\text{mm}$



SLA





PCB



Large inductor
for RF drive

Small amplifier
for RF drive

Arrays offset the effect of ion attenuation due to miniaturization



- Fabrication Technology:
Stereolithography apparatus
- multiple RITs for increased ion capacity
- Geometrically complex electrode structures
 - Circular
 - Linear
- Fabrication flexibility
 - Monolithic construction
 - Independently isolated electrodes




Stereolithography Apparatus (SLA) Fabrication

Trapping dimensions:

Analyzer footprint:

Base PCB substrate:

Viper si2 SLA system






Additive Process
Monolithic, fully integrated
analyzer


Flexible placement of
electrodes & connections

Laser ablation used for
electrode isolation

Array Operation and Control







Array Operation and Control: Bias Tee as a On & Off Switch





- Selective ion trap(s) operation inside the array
- Fast, reversible, and outside the vacuum chamber
- No physical obstruction to ion injection or detection
- RF voltages applied to all traps in the array


Performance: Full scan MS (SLA vs. SS RIT)






Performance: Ion injection control, best resolution, and Intact protein analysis





Performance: Tandem MS of atenolol, double DAPI, Full, Isolation, and MS2



isolation q: 0.434
SWIFT notch: 148.7 to 154.9
isolation time: 75 ms
ion cooling time: 1750 ms
activation q: 0.655
activation time: 30 ms
activation amp: ~1 V_{p-p}
DAPI time: 30 ms
KVDAC: 1.68 Vdc
dDAPI time: 10 ms




Multi-channel data: All eight channels





- RF voltages applied to all traps during analysis
- Primary cause of mass shift due Δr_0
- Ejection occurs at a $\Delta \sim 5 V_{0-p}$ for same m/z
- Δr_0 of $\sim 10\text{s}$ of μm between adjacent traps
- Compensate using individual RF drives

Multi-channel data: Response for multiple vs. single channel





- 0.1 $\mu\text{m}/\text{ml}$ to 10 $\mu\text{g}/\text{ml}$ dynamic range
- Increased sensitivity & dynamic range with multiple channels
- Development of custom peak detection algorithm
- Peak widths for combined channels 2 – 2.5 amu without RF compensation


Multi-channel RF drive: High Performance Operational Amplifiers




- Independent RF for each ion trap in the array
- Apex PA85 operational amplifier
- 450 V_{p-p} output maximum voltage
- 1000 V/ μ s slew rate
- Anticipate maximum m/z 500 at 950 kHz
- Cost effective to produce 4 layer PCB




Multi-channel RF drive: Fun with Independent RF Control



- *N,N-Diethyl-3-methylbenzamide*
- Scan trapping 1 before 4 of channel 1 relative to 4



Combined technology: Ion trap array design, fabrication, and control electronics



Acknowledgements



Aston Laboratory

Abraham Badu-Tawiah
Brandy Houmes
Dahlia Campbell
Fatkhulla
Tadjimukhamedov
Graham Cooks
Guangming Huang
Guangtao Li
He Wang
Isabella Zhang

Jiangjiang Liu
Jobin Cyriac
Jon Dalgleish
Josh Wiley
Leann Herndon
Livia Eberlin
Nick Manicke
Ryan Espy
Sandilya Garimella
Santosh Soparawalla

Sheran Oradu
Kevin Kerian

Center for Analytical Instrumentation Development



Matt McNickels
Jason Duncan
Matt Kirleis
Frank Boudreau
Robert Noll

Thermo Fisher Scientific:
132003660

DHS grant:
DHSHSHQDC-09-9-00008

NSF grant:
CHE0848650

Questions

SLA Fabrication: Build Process and Laser Ablation



Ion Trap Build Process



Laser Ablation: Electrode Isolation